

Von Kossa Stain

7/27/2023

Slide with histological specimen	
Distilled water	1 min
Reagent 1 (Silver nitrate solution) under exposure to illuminant	120 min
Running tap water	3 min
Reagent 2 (Sodium thiosulfate solution)	5 min
Running tap water	1 min
Nuclear fast red-aluminium sulfate solution 0.1% (optional)	7 sec
Distilled water	1 min
Ethanol 70 %	1 min
Ethanol 96 %	1 min
Ethanol 100 %	1 min
Ethanol 100 %	1 min
Xylene	5 min
Xylene	5 min
Mount the slides with Perimount and cover glass before Xylene dries.	